



Custom Ear Plugs

User Instructions



Please note that all our custom hearing products come with a 30-Day Fit & Filter Guarantee to ensure you are completely satisfied with your purchase.



For more information contact:
Universal Aids Ltd

8-14 Wellington Road South, Stockport, UK, SK4 1AA

Tel: +44 (0) 161 480 9228 Fax: +44 (0) 161 476 5707

Email: protection@db-hearing.co.uk

Web: www.db-hearing.co.uk

SnorePlug

CUSTOM HEARING PROTECTION

SnorePlug

Re-Usable Earplugs

User Instructions

Fitting Instructions

Insertion: Firstly ensure that you have the correct side earplug for the ear. This can normally be identified by the coloured text on the outside of the plug (Blue—Left & Red—Right). The plugs may also be coloured differently for each ear with the colour choice having been agreed with your audiologist.

Check that the earplug is clean then place the earplug text-side-down on a flat surface or the palm of your hand (Fig 1), hold the plug between the thumb and forefinger and apply a small amount of the gel (supplied) on to the ear tip (Fig 2). Now holding the plug firmly, place the canal tip into the ear canal (Fig 3). The ear plug can then be firmly seated into place by applying gentle pressure with the little finger. When fitted correctly the earplug will fill the recesses of the ear and feel snug and the removal tab will be lying flush in the bowl of the ear (Fig 4).

In order to remove the plug, grip the removal tab using finger and thumb in the bowl of the ear and gently pull forward to break the seal. The plug can then be removed, wiped clean and returned to the case.

Take care not to apply excessive force to the Tab Handle.

Cleaning and Disinfecting: The earplug should be wiped over with a soft damp wipe or cloth and dried after use. SnorePlug can be disinfected using a very mild solution such as Milton™.

Storage: SnorePlug comes with a neat pocket size Glow-in-the-Dark case which is ideal for storing your sleep plugs. The earplugs should be inspected regularly for damage and always stored in the case before and after use.

IMPORTANT: SnorePlug is designed to be a comfort plug which helps to reduce sound such as snoring, enabling the user to have a relaxing sleep. It is not intended for use in extremely high levels of noise where suitable hearing protection should be worn. SnorePlug aims to combine the reduction of sound alongside the comfort of a custom fit design, and there is therefore no guarantee that all sound will be blocked out while wearing this product. Care should be taken not to use excessive force on the removal tab or use fingernails to grip the item.

