




Earmould Tubing and Accessories




Standard Wall Pre Bent Earmould Tubing made from clear, medical grade PVC. Outside diameter 3mm, inside diameter 2mm. Available in bags of 100
Universal dB Product Code: 221
NHS Code: GHD099



Double Bend Pre Bent Earmould Tubing made from clear medical grade PVC. Outside diameter 3mm, inside diameter 2mm. Available in bags of 100
Universal dB Product Code: 349
NHS Code: GHD2735



Thick Wall Pre Bent Earmould Tubing made from clear medical grade PVC. Outside diameter 3.3mm, inside diameter 2mm. Available in bags of 100
Universal dB Product Code: 222



Low Condensate-Stay dry Pre Bent Earmould tubing which is designed to reduce moisture build up due to condensation for hearing aid wearers
Universal dB Product Code: 263



Tube Lock Tubing with barbed collar which helps retain tubing within soft silicone type earmoulds.
Universal dB Product Code: 264



Libby Horn Tubing (3mm or 4mm)
Universal dB Product Code: 152



Coloured Pre Bent Tubing (Standard Wall) available in the following colour options:
Red, Blue, Yellow, Green, Black & Brown
Universal dB Product Code: 211



Standard Wall /Thick Wall Earmould Tubing made from clear, medical grade PVC. (100m Roll)
Universal dB Product Codes:
Standard Wall Code 215
Thick Wall Code 214



Barbed Elbow Tubing Adaptor for earmoulds
Universal dB Product Code: 240



Tube Lock Insert retention aid for tubing, for use in all soft earmoulds (gold plated or plastic options)
Universal dB Product Code: 264



Select-a-Vent tree. Range of vent options in various sizes to fit earmould vents.
Universal dB Product Code: 272



Cedis Skin Care Gel with Aloe Vera used to prevent irritation and ease of insertion of earmoulds, swim moulds and hearing protection
Universal dB Product Code: 347



Audio Spa earmould hygiene system with Nitrodine™ soluble sterilizing tablets.
Highly effective at removing microbial biofilm from synthetic polymers
Universal dB Product Code: 292