

CODE

EM4112
EM4121
EM4112 Semi
EM4121 Semi
EM2108 BIO ST
EM2108 BIO HT
EM2109 BIO ST
EM2109 BIO HT

DIGITAL BTE'S – Hypoallergenic

Ultra Carved Shell, Medium Hardness Biopore Digital **HP**
Ultra Carved Skeleton, Medium Hardness Biopore Digital **HP**
Semi Carved Shell, Medium Hardness Biopore Digital **HP**
Semi Carved Skeleton, Medium Hardness Biopore Digital **HP**
Biopore shell with soft tip **HP**
Biopore shell with hard tip **HP**
Biopore skeleton with soft tip **HP**
Biopore skeleton with hard tip **HP**

Biopore is also available in soft & hard shore hardness
(clear colour only) by request.

'CLICK' - CUSTOM OPEN Moulds (For Corda Type Micro Tube)

Please enquire for more details

EM2800 Custom Open CLICK - Hard Tip
EM2800 CL Custom Open CLICK - Hard Tip with canal lock
EM2900 Custom Open - Soft Tip
EM2900 CL Custom Open Fitting - Soft Tip with canal lock

EARPHONE ADAPTORS

EM2180 PreciseSound Custom Earphone earmoulds (many options available)

EARMOULD OPTIONS

218 Thick Wall Pre Bent Tube
290 Double Bend Pre Bent Tube
247 Baby Wire
287 Removal Handle
240 Barbed Elbow
272 Select-A-Vent
211 Coloured Pre Bent Tube
150 High Power
238 Glitter
239 Picture Insert
263 Low Condensate Tube
262 Ultra Carved Profile
261 Marble Effect
153 Acoustic Filter
152 Libby Horn
282 Stay Dry Libby Horn
264 Tube Lock Tubing
275 Non Allergic Coating
250 Swim Plug Cord
177 Gold Plating
324 Nanoscreen™ Anti Bacterial Coating (for hard and silicone earmoulds)

Universal dB

CUSTOM HEARING PRODUCTS



EARMOULD REFERENCE GUIDE

Universal Aids Ltd
8-14 Wellington Road South
Stockport, Cheshire,
England SK4 1AA

T: +44 (0)161 480 9228

F: +44 (0)161 476 5707

E: hearing@universaldb.co.uk

CODESOLID (Body Worn) Standard

EM2101	Hard Acrylic
EM2103	Hard Acrylic with Baby Wire
EM2104	Hard Acrylic Soft Tip
EM2105	Hard Acrylic Soft Tip with Baby Wire
EM2113	Soft Acrylic with Hard Back
EM2120	Soft Acrylic with Hard Back and Baby Wire

SOLID (Body Worn) – Non Allergic Materials

EM2102 – ACR	Acrylon - Urethane alternative to acrylic
EM2102 – AVC	Soft Acrylic with hard Softens at body temperature for an superior fit. HP*
EM2102 – Bio	Biopore - Silicone with gold plated clip(available in standard, hard or soft)
EM2102 – MF	Microflex - Soft silicone with gold plated clip
EM2102 – NC	BioMed Non Cadmium –Digital Hard acrylic
EM2102 – VSA	Soft acrylic

SHELL Standard

EM2106	Hard acrylic
EM2110	Hard acrylic soft tip
EM2112	Soft acrylic

SHELL – Non Allergic Materials

EM2108 – ACR	Acrylon - Urethane alternative to hard acrylic
EM2108 – AVC	Soft acrylic , softens at body temperature for an superior fit. HP
EM2108 – Bio	Biopore silicone (available in standard, hard or soft) HP
EM2108 – MF	Microflex - Soft silicone HP
EM2108 – NC	BioMed Non Cadmium –Digital Hard UVcure acrylic
EM2108 – VSA	Soft Acrylic , softens at body temperature for an improved fit.
EM2108 – THERM	ThermoTEC Thermal flexible Polyurethane material HP

NANO Earmould Coating

NANO 324	NANO Anti-Bacterial coating is available for all hard earmould materials in our range. A NANO coating is also available for soft, silicone type earmoulds including Biopore and Microflex. Unfortunately NANO is unavailable for soft acrylic moulds such as 2112, 2108 VSA and soft tip 2111 &2110 earmoulds.
----------	--

***HP** Indicates suitable for use with High Power hearing aids

CODESKELETON Standard

EM2107	Hard Acrylic
EM2111	Hard Acrylic Soft Tip
EM2117	Hard Acrylic Half Phantom
EM2118	Hard Acrylic Soft Tip Half Phantom
EM2121	Soft Acrylic
EM2122	Open Mould - Retention By Ring, Claw or Gap (see open mould chart)

SKELETON Non Allergic Materials

EM2109 - ACR	Acrylon - Urethane alternative to hard acrylic
EM2109 - AVC	Soft Acrylic , softens at body temperature for an superior fit.
EM2109 – Bio	Biopore silicone Digital (available in standard, hard or soft)
EM2109 – MF	Microflex - Soft silicone
EM2109 – NC	BioMed Non Cadmium –Digital Hard UV cure acrylic
EM2109 – S	Soft Acrylic
EM2109 – VSA	Soft Acrylic , softens at body temperature for an improved fit.
EM2109 – THERM	ThermoTEC Thermal flexible Polyurethane material
EM2122 – THERM	ThermoTEC Thermal flexible Polyurethane material

CANAL Earmoulds Standard Earmould Tube

EM2114	Meatal Tip, hard acrylic
EM2115	Meatal Tip, hard acrylic with soft tip
EM2116	Meatal Tip, soft acrylic
EM2114 CL	Meatal Tip, hard acrylic with canal lock
EM2115 CL	Meatal Tip, hard acrylic with soft tip with canal lock
EM2116 CL	Meatal Tip, soft acrylic with canal lock

WATER PROTECTION

EM2119	Swim Plug soft acrylic
EM2119 – CAT	Swim Plug, floatable
EM2119 – MF	Swim Plug, Microflex
EM3119	WaterBarrier™ - Floatable, Anti Bacterial Swim Moulds With Case , Skin Care Gel & Care Instructions
EM5119	WaterBarrier™ SURF WaterBarrier designed for Watersports

HEARING PROTECTION

EM2365	SnorePlug™- Soft Canal Ear Plug with Case and Skin Care Gel
EM2390	NoiseBarrier™ - Certified Custom Hearing Protection with Filter
EM2321	ER Series 9/15/25 Noise Protector with Filter, for musicians
EM2321BANG	ER BANG Filtered protection for Shooters
EM2395	CENS Proflex hearing protection for shooting
EM2360	Unfiltered Noise Plug, Biopore
EM9300/9400/9500	MiniEar Custom protection for Musicians/Industry/shooting