

## CODE

EM4112  
EM4121  
EM4112 Semi  
EM4121 Semi  
EM2108 BIO ST  
EM2108 BIO HT  
EM2109 BIO ST  
EM2109 BIO HT

## DIGITAL BTE'S – Hypoallergenic

Ultra Carved Shell, Medium Hardness Biopore **Digital HP**  
Ultra Carved Skeleton, Medium Hardness Biopore **Digital HP**  
Semi Carved Shell, Medium Hardness Biopore **Digital HP**  
Semi Carved Skeleton, Medium Hardness Biopore **Digital HP**  
Biopore shell with soft tip **Digital HP**  
Biopore shell with hard tip **Digital HP**  
Biopore skeleton with soft tip **Digital HP**  
Biopore skeleton with hard tip **Digital HP**

Biopore is also available in soft & hard shore hardness  
(clear colour only) by request.

## 'CLICK' - CUSTOM OPEN Moulds (For Corda Type Micro Tube)

Please enquire for more details

EM2800 Custom Open CLICK - Hard Tip **Digital**  
EM2800 CL Custom Open CLICK - Hard Tip with canal lock **Digital**  
EM2900 Custom Open - Soft  
EM2900 CL Custom Open Fitting - Soft with canal lock

## EARPHONE ADAPTORS

EM2180 PreciseSound Custom Earphone earmoulds (many options available)

## EARMOULD OPTIONS

218 Thick Wall Pre Bent Tube  
290 Double Bend Pre Bent Tube  
247 Baby Wire  
287 Removal Handle  
240 Barbed Elbow  
272 Select-A-Vent  
211 Coloured Pre Bent Tube  
150 High Power  
238 Glitter  
239 Picture Insert  
263 Low Condensate Tube  
262 Ultra Carved Profile  
261 Marble Effect  
153 Acoustic Filter  
152 Libby Horn  
282 Stay Dry Libby Horn  
264 Tube Lock Tubing  
275 Non Allergic Coating  
250 Swim Plug Cord  
177 Gold Plating  
324 Nanoscreen™ Anti Bacterial Coating (for hard, and silicone earmoulds)

# Universal dB

CUSTOM HEARING PRODUCTS



## EARMOULD REFERENCE GUIDE

Universal Aids Ltd  
8-14 Wellington Road South  
Stockport, Cheshire,  
England SK4 1AA

T: +44 (0)161 480 9228  
F: +44 (0)161 476 5707  
E: [hearing@universaldb.co.uk](mailto:hearing@universaldb.co.uk)

E: [hearing@universaldb.co.uk](mailto:hearing@universaldb.co.uk) W: [www.db-hearing.co.uk](http://www.db-hearing.co.uk)

CODESOLID (Body Worn) Standard

EM2101	Hard Acrylic
EM2103	Hard Acrylic with Baby Wire
EM2104	Hard Acrylic Soft Tip
EM2105	Hard Acrylic Soft Tip with Baby Wire
EM2113	Soft Acrylic with Hard Back
EM2120	Soft Acrylic with Hard Back and Baby Wire

SOLID (Body Worn) – Non Allergic Materials

EM2102 – ACR	Acrylon - Urethane alternative to acrylic
EM2102 – AVC	Soft Acrylic Softens at body temperature for an superior fit. <b>HP*</b>
EM2102 – Bio	Biopore - Silicone with gold plated clip standard, hard or soft shore)
EM2102 – MF	Microflex - Soft silicone with gold plated clip
EM2102 – NC	BioMed Non Cadmium – <b>Digital</b> Hard acrylic
EM2102 – VSA *	Soft acrylic

SHELL Standard

EM2106	Hard acrylic
EM2110	Hard acrylic soft tip
EM2112	Soft acrylic

SHELL – Non Allergic Materials

EM2108 – ACR	Acrylon - Urethane alternative to hard acrylic
EM2108 – AVC	Soft acrylic , softens at body temperature for an superior fit. <b>HP</b>
EM2108 – Bio	Biopore silicone <b>Digital</b> (available in standard, hard or soft shore) <b>HP</b>
EM2108 – MF	Microflex - Soft silicone <b>HP</b>
EM2108 – NC	BioMed Non Cadmium – <b>Digital</b> Hard UVcure acrylic
EM2108 – VSA *	Soft Acrylic , softens at body temperature for an improved fit.
EM2108 – THERM	ThermoTEC Thermal flexible Polyurethane material <b>HP</b>

NANO Earmould Coating

NANO 324	NANO Anti-Bacterial coating is available for all hard earmould materials in our range. A NANO coating is also available for soft, silicone type earmoulds including Biopore and Microflex. Unfortunately NANO is unavailable for soft acrylic moulds such as 2112, 2108 VSA and soft tip 2111 & 2110 earmoulds.
----------	---

\***HP** Indicates suitable for use with High Power hearing aids

\* Please note that if a Non Allergic material is not specified the the default VSA material will be used.

E: hearing@universaldb.co.uk W: www.db-hearing.co.uk

CODESKELETON Standard

EM2107	Hard Acrylic
EM2111	Hard Acrylic Soft Tip
EM2117	Hard Acrylic Half Phantom
EM2118	Hard Acrylic Soft Tip Half Phantom
EM2121	Soft Acrylic
EM2122	Open Mould - Retention By Ring, Claw or Gap (see open mould chart)

SKELETON Non Allergic Materials

EM2109 - ACR	Acrylon - Urethane alternative to hard acrylic
EM2109 - AVC	Soft Acrylic , softens at body temperature for an superior fit.
EM2109 – Bio	Biopore silicone <b>Digital</b> (available in standard, hard or soft shore)
EM2109 – MF	Microflex - Soft silicone
EM2109 – NC	BioMed Non Cadmium – <b>Digital</b> Hard UV cure acrylic
EM2109 – S	Soft Acrylic
EM2109 – VSA *	Soft Acrylic , softens at body temperature for an improved fit.
EM2109 – THERM	ThermoTEC Thermal flexible Polyurethane material
EM2122 – THERM	ThermoTEC Thermal flexible Polyurethane material

CANAL Earmoulds Standard Earmould Tube

EM2114	Meatal Tip, hard acrylic
EM2115	Meatal Tip, hard acrylic with soft tip
EM2116	Meatal Tip, soft acrylic
EM2114 CL	Meatal Tip, hard acrylic with canal lock
EM2115 CL	Meatal Tip, hard acrylic with soft tip with canal lock
EM2116 CL	Meatal Tip, soft acrylic with canal lock

WATER PROTECTION

EM2119	Swim Plug soft acrylic
EM2119 – CAT	Swim Plug, floatable <b>Digital</b>
EM2119 – MF	Swim Plug, Microflex
EM3119	WaterBarrier™ - Floatable, Anti Bacterial Swim Moulds With Case , Skin Care Gel & Care Instructions <b>Digital</b>
EM4119	WaterBarrier™ <b>Alert!</b> Water protection Featuring <b>Acoustic Filter Digital</b>
EM5119	WaterBarrier™ <b>SURF</b> WaterBarrier designed for Watersports <b>Digital</b>

HEARING PROTECTION

EM2360	Unfiltered Noise Plug, Biopore <b>Digital</b>
EM2365	SnorePlug™- Soft Canal Ear Plug with Case and Skin Care Gel <b>Digital</b>
EM2390	NoiseBarrier™- Certified Custom Hearing Protection with Filter <b>Digital</b>
EM2321	ER Series 9/15/25 Noise Protector with Filter, for musicians <b>Digital</b>
EM9300/9400/9500	MiniEar™ Custom protection for Music/Industry/Shooting <b>Digital</b>

E: hearing@universaldb.co.uk W: www.db-hearing.co.uk